

Colony PCR

1. Pick a colony, make a replica on LB-plate with antibiotic. Transfer bacteria to 50 μ l of sterile water into 0.5ml tube. Vortex.
2. Boil for 5min. Vortex again.
3. Spin @ 11,000g for 1min to pellet down cell debris.
4. Use 10 μ l of supernatant for PCR (keep on ice).
5. Make a master mix for PCR:

MASTERMIX (multiply by number of analyzed colonies + a negative control: a colony from negative control plate):

H ₂ O	7.05
10xTAQ buffer	2.5
DMSO (5% final)	1.25
25mM MgCl ₂	1.0
25mM dNTPs	0.2
10 μ M up primer	1.25
10 μ M dn primer	1.25
TAQ polymerase	<u>0.5</u>
	15 μ l

Do NOT use high fidelity polymerase (such as Phusion)!

6. Add 15 μ l of PCR master mix to each tube containing 10 μ l of colony supernatant.

Program: COLONY

- 95°C – 2min
- 35 cycles:
 - 95°C – 30sec
 - primer-specific anneal. T°C – 30sec
 - 72°C – (1min per kb)
- 72°C – 10min
- 4°C.