

Site Directed Mutagenesis Protocol

Based on the QuikChange™ Site Directed Mutagenesis Kit

1. Follow the QuikChange™ recommendations for primer design:
 - a. Both mutagenic primers must contain the desired mutation and anneal to the same sequence on opposite strands of the plasmid
 - b. Primers should be 25-45 bases in length
 - c. T_m should be greater than or equal to 78°C

The following formula is commonly used for estimating the T_m of primers:
 $T_m = 81.5 + 0.41(\%GC) - 675/N - \% \text{ mismatch}$, where

 - N is the primer length in bases
 - values for % GC and % mismatch are whole numbers

For calculating T_m for primers intended to introduce insertions or deletions, use this modified version of the above formula:
 $T_m = 81.5 + 0.41(\%GC) - 675/N$, where N does not include the bases, which are being inserted or deleted.
 - d. Desired mutation must be in the middle of the primer with ≥ 15 bp on both sides of the mutation/insertion/deletion
 - e. CG content should be higher than 40%
 - f. 3'-end of primers should end in 1 or more G or C
 - g. Protocol says they must be HPLC purified but when ordering, desalted have always worked just fine
2. For each 50 μ l reaction, use the following protocol for the reactions and add in the order listed below:

Stock	Final	Content	Volume (μ l)
		ddH ₂ O	Bring to 50 μ l
Tube at -20	1 M	Ethylene Glycol	2.8
	25 ng	Template DNA	
5X	1X	Phusion Buffer	10
25 mM	0.4 mM	dNTPs	0.8
20 μ M	0.2 μ M	Forward Primer	0.5
20 μ M	0.2 μ M	Reverse Primer	0.5
		*Phusion Polymerase	0.5

*Polymerase should not be added to the reaction until the thermocycler has begun heating up to 98C

3. Program the Thermocycler to the following set up for the reaction:

Cycle	Temperature	Time	Repeats
1	98°C	10 sec	1
2	98°C	30 sec	12-18 Times [#]
	55°C	60 sec	
	72°C	30 sec/kb	
3	4°C	∞	

#	Type of mutation desired	Number of cycles
	Point mutations	12
	Single amino acid changes	16
	Multiple amino acid deletions or insertions	18

4. Add 1 μ l of DpnI endonuclease to the PCR tube and place at 37°C for at least 2 hours.
5. Transformation.